



# Large **V-Wipes** Alcohol Free Surface Disinfectant

For the surface disinfection of non-invasive medical devices, including surfaces sensitive to alcohol

## Key Product Features

Effective against a wide spectrum of organisms (see bullets opposite). Tubs contain 200 wipes in polysealed bags. Wipe size 200 x 200mm. Dentomint fragrance. Also contains detergent.

## Directions For Use

Using a new wipe for each application, wipe objects to be disinfected so that they are completely moistened and let the disinfectant take action over a period of 30 seconds, allowing surface moisture to evaporate. Wipe with clean tissue if necessary.

## Product Composition

Solution contains: Aqua, Benzalkonium Chloride, C12-14 Alkyl Ethylbenzylammonium Chloride, Didecylmethylammonium Chloride, Isopropyl Alcohol, C9-11 Pareth-8, Peracetic acid, PPG-2 Methyl Ether,

Polyaminopropyl-Biguanide, Cocamidopropyl Betaine, Tetrasodium EDTA.

## Physical Properties

PolyViscose Substrate impregnated with clear liquid. Colour white. Odour mint. pH 8-9 of solution with a shelf life of 3 years.

## Test Standards

Tested to standards  
BS EN14476  
BS EN13727  
BS EN14204  
BS EN13704  
BS EN13624

Also tested to DGHM guidelines and AFNOR standards.

## Hazard Classifications and Details

This product is not classified as dangerous for carriage.

## Packing Details

12 x Tub of 200 wipes in a box  
Size 340 x 450 x 255 (h) mm  
Weight per box 8.4kg  
12 x Refill of 200 wipes in a box  
Size 250 x 300 x 295 (h) mm  
Weight per box 7.2kg

## Packaging and Accessories

BOS9810 Tub of 200 large wipes  
BOS9811 Wipe refill pack of 200 large wipes

## Effective Spectrum

### Bacteria

- MRSA
- Escherichia coli
- Pseudomonas aeruginosa
- Enterococcus hirae
- Staphylococcus aureus
- Salmonella typhimurium
- Listeria monocytogenes
- Legionella pneumophila
- Serratia marcescens
- Vancomycin-resistant enterococcus faecium
- Proteus mirabilis
- Klebsiella pneumoniae
- Mycobacterium tuberculosis Avium
- Mycobacterium tuberculosis Terrae
- Mycobacterium bovis

### Viruses

- Hepatitis B virus
- Hepatitis C virus
- HIV (HIV- type 1)
- Herpes simplex virus type 1
- H1N1 Influenza A virus
- Rotavirus A

### Fungus

- Aspergillus niger
- Candida albicans
- Penicillium chrysogenum
- Penicillium verrucosum

### Spores

- Clostridium difficile (C-Diff)
- Cladosporium

CE  
0120



“ 200 x 200mm  
Large wipes ”